

Sensitive detection of circulating tumor DNA by whole genome sequencing

Validation of MRDetect using serial blood samples from stage III colorectal cancer patients

Amanda Frydendahl¹, Jesper Nors¹, Thomas Reinert¹, Sunil Deochand², Dillon Maloney², Noah Friedman², Tomer Lauterman³, Danielle Afterman³, Imane Bourzgui², Nidhi Ramaraj², Zohar Donenhirsh³, Ronel Veksler³, Ravi Kandasamy², Iman Tavassoly², Jonathan Rosenfeld², Anders Husted Madsen⁴, Uffe S. Løve⁵, Per Vadgaard Andersen⁶, Ole Thorlacius-Ussing⁷, Lene Hjerrild Iversen⁸, Kåre Andersson Gotschalck⁹, Boris Oklander³, Asaf Zviran², Claus Lindbjerg Andersen¹

¹Department of Molecular Medicine, Aarhus University Hospital, Denmark, ²C2i Genomics Inc., New York, USA, ³C2i Genomics, Ltd., Haifa, Israel, ⁴Department of Surgery, Herning Regional Hospital, Denmark, ⁵Department of Surgery, Viborg Regional Hospital, Denmark, ⁶Department of Surgery, Odense University Hospital, Denmark, ⁷Department of Surgical Gastroenterology, Aalborg University Hospital, Denmark, ⁸Department of Surgery, Aarhus University Hospital, Denmark, ⁹Department of Surgery, Randers Regional Hospital, Denmark



Amanda Frydendahl Johansen

amandafbj@clin.au.dk

LinkedIn:

Conclusions

- 1) Detection of circulating tumor DNA (ctDNA) through whole genome sequencing (WGS) of plasma is associated with high risk of recurrence
- 2) Our approach show great reproducibility across independent laboratories
- 3) Detection of ctDNA immediately after surgery may be affected by high amounts of trauma-induced circulating free DNA (cfDNA)
- 4) Differences between tumor and plasma WGS suggest tumor heterogeneity

Study summary

Background

- Detection of circulating tumor DNA (ctDNA) is associated with poor prognosis.

- Sensitive detection of ctDNA is challenged by very low amounts of ctDNA, often only one mutated copy in 10 mL blood.

- We developed MRDetect; a whole genome sequencing approach, which detects ctDNA
- using a patient-specific signature generated from tens of thousands of genetic alternations throughout the genome.

Aim


- Validation of MRDetect for sensitive detection of ctDNA to monitor disease in stage III colorectal cancer patient, including assessment of reproducibility.

Methods

- Reproducibility: Analysis of paired aliquots of 136 blood samples from 14 patients.

- Clinical utility of MRDetect: Whole genome sequencing (WGS) of tumor (n = 129), normal (n = 129) and serially collected plasma (n = 921) samples for ctDNA assessment using MRDetect.


Workflow




Results

Prognostic impact of post-OP ctDNA status

Timing of post-OP sampling affects detection of ctDNA




High cfDNA levels post-OP increase the risk of false negative ctDNA detection



Disease surveillance using ctDNA

Serial assesment of ctDNA idenfies early recurrence

A total of 102 patients (16 recurrence, 86 non-recurrence) are included in the surveillance analysis



WGS of tumor and plasma reveals tumor heterogeneity

Chromosome 3 copy number change adjusting mutant SETD2 frequency






Figure 7
Case study indicating tumor heterogeneity detected through whole genome sequencing (WGS) of plasma samples. WGS of metastasis confirms findings in plasma.

Study population



Reproducibility



Post-ACT ctDNA status and prognosis

ctDNA assesment after adjuvant chemotherapy identifies patients with residual disease

