

Clinical Validation of the Afirma Genomic Sequencing Parathyroid Classifier

Julie Ann Sosa,¹ Trevor E. Angell,² Joshua Barbiariz,³ Neil Barth,³ Thomas Blevins,⁴ Quan-Yang Duh,⁵ Ronald A. Ghossein,⁶ R. M. Harrell,⁷ Jing Huang,³ Urooj Imtiaz,³ Giulia Kennedy,³ Su Yeon Kim,³ Richard T. Kloos,³ Virginia A. LiVolisi,⁸ Keping N. Patel,⁹ Gregory Randolph,¹⁰ Peter M. Sadow,¹¹ Michael H. Shanik,¹² S. T. Traweck,¹³ P. S. Walsh,³ Duncan Whitney,³ Michael Yeh,¹⁴ Paul W. Ladenson¹⁵


(1) Department of Surgery, Section of Endocrine Surgery, Duke University Medical Center, Durham, NC (2) Department of Medicine, Division of Endocrinology, Diabetes, and Hypertension, Brigham and Women's Hospital and Harvard Medical School, Boston, MA (3) Veracyte, Inc., South San Francisco, CA (4) Texas Diabetes and Endocrinology, Austin, TX (5) Department of Surgery, Section of Endocrine Surgery, University of California San Francisco, San Francisco, CA (6) Department of Pathology, Division of Head and Neck Pathology, Memorial Sloan-Kettering Cancer Center, New York, NY (7) The Memorial Center for Integrative Endocrine Surgery, Boca Raton, FL (8) Department of Pathology and Laboratory Medicine, Anatomic Pathology Division, University of Pennsylvania School of Medicine, Philadelphia, PA (9) Department of Surgery, Division of Endocrine Surgery, NYU Langone Medical Center, New York, NY (10) Department of Otolaryngology, Division of Thyroid and Parathyroid Endocrine Surgery, Massachusetts Eye and Ear and Harvard Medical School, Boston, MA (11) Department of Pathology, Massachusetts General Hospital and Harvard Medical School, Boston, MA (12) Endocrine Associates of Long Island, Smithtown, NY (13) Thyroid Cytopathology Partners, Austin, TX (14) Department of Surgery, Endocrine Surgery Program, UCLA David Geffen School of Medicine, Los Angeles, CA (15) Department of Medicine; Division of Endocrinology, Diabetes and Metabolism, Johns Hopkins University, Baltimore, MD



INTRODUCTION

The parathyroid glands are located adjacent to the thyroid and occasionally within it. Enlarged and imbedded parathyroid glands can be mistaken as thyroid nodules or suspicious lymph nodes. On fine needle aspiration biopsy (FNAB) of such lesions, cytology is often indeterminate, failing to identify its parathyroid origin and potentially resulting in an unnecessary thyroid surgery. The Afirma Genomic Sequencing Classifier (GSC) identifies genomically benign thyroid nodules among those with indeterminate FNAB to prevent unnecessary diagnostic surgery using RNA sequencing and machine learning algorithms. Integrated cassettes are used to detect the molecular signatures of specific neoplasms that may further alter patient care, such as unsuspected parathyroid tissue (Figure 1). Here we report the clinical performance of the parathyroid classifier cassette used with the GSC.

FIGURE 1. Afirma Genomic Sequencing Classifier (GSC) Integrated Workflow



METHODS

Algorithm training was performed with a set of 476 FNAs—6 parathyroid and 470 thyroid FNAs. An additional 97 tissues were used in feature selection, but not model training (Table 1). A support vector machine classifier was developed using 109 differentially expressed genes, including the 5 genes in the Afirma GEC parathyroid cassette (Figure 2).

TABLE 1. TRAINING COHORTS

1A. Core FNA Training Set (Parathyroid=6, Non-Parathyroid=470)

Parathyroid Cohort

Cohort	Bethesda III	Bethesda IV
ENHANCE Arm1	2	0
Parathyroid	3	1

Non-Parathyroid Cohort


Cohort	Bethesda III	Bethesda IV
ENHANCE Arm1	209	76
ENHANCE Arm2	50	14
Cytol GECB	111	
MTC	5	5

1B. Auxiliary Tissue Training Set (n=97)

Cohort	Parathyroid	Non-Parathyroid
MTC	0	21
Parathyroid	21	0
Hürthle Tissue	0	25
Tissue	0	30

FIGURE 2. Relative Importance of 109 Parathyroid Classifier Genes

109 genes used in the classifier are sorted by their relative importance where 100% indicates the highest importance. 5 genes that were used in the previous version of the parathyroid cassette are highlighted in red in this figure.



■ 5 genes in the Afirma GEC parathyroid cassette
■ The other 104 differentially expressed genes used in the classifier

RESULTS

The final locked classifier was blindly validated on an independent test set of 195 FNAs (118 Bethesda III, 77 Bethesda IV). The classifier had 100% sensitivity [4/4 parathyroid correctly called positive; CI 39.8-100%] and 100% specificity [191/191 thyroid correctly called negative; CI 98.1-100%] (Table 2). All positive samples had clinical/surgical confirmation of the parathyroid etiology, while all negative samples were negative on surgical pathology.

TABLE 2. Independent Validation Performance (n=195)

		Truth Label (n = 195)	
		Parathyroid	Non-Parathyroid
Parathyroid Classifier Result	Positive	4	0
	Negative	0	191
		Sensitivity	Specificity
		100% (39.8-100%)	100% (98.1-100%)

CONCLUSIONS

Preoperative genomic identification of parathyroid tissue facilitates appropriate management of parathyroid tissue within and adjacent to the thyroid gland. Reporting the result of the Parathyroid Classifier with every GSC test may help avoid inappropriate interventions, costs, and potential complications.¹

1. Meltzer, C et al. *Otolaryngology–Head and Neck Surgery* 2016.