The Genomic Landscape of Thyroid Nodules in Patients ≥75 Years Old





Architha Kannan, MBBS,¹ Yang Chen, PhD,² Mohammed Alshalalfa, PhD,² Joshua P. Klopper, MD,² Sanziana A. Roman, MD¹

1. Division of Surgical Oncology, Section of Endocrine Surgery, University of California San Francisco School of Medicine, San Francisco, CA. 2. Veracyte, Inc., South San Francisco, CA.

INTRODUCTION

- Thyroid nodules are increasingly common in older adults; most are benign.¹
- In patients ≥75 years, malignant nodules often behave more aggressively, making accurate risk stratification essential.
- The Afirma Genomic Sequencing Classifier (GSC) is an exome-enriched RNA-seq assay validated in adults ≥21 years that refines malignancy risk in Bethesda III/IV nodules, potentially avoiding unnecessary surgery.²
- Limited data exist on the genomic and transcriptomic landscape of thyroid nodules in patients ≥75 years.
- The goal of this study was to characterize the molecular and transcriptomic profiles of thyroid nodules in patients
 ≥75 years compared with younger cohorts.

METHODS

- The Afirma GSC database was queried from August 2017 June 2024.
- Nodules that had adequate genomic material from follicular cells were identified to analyze expressed variants and fusions among Afirma GSC-(S)uspicious and nodules with Bethesda V/VI cytology.
- Nodules from patients <20 yrs were excluded.
- TERT promoter mutations were analyzed from March 2023 – June 2024.
- A subset of samples with available transcriptomic were analyzed for ERK and EMT signaling expression, the BRAF-RAS (BRS) score, and expression of the Human Molecular Signatures Database (MSigDB) hallmarks of cancer pathways.
- Patients were stratified into the following age cohorts:
- 20-54 yrs
- 55-74 yrs
- ≥75 yrs

RESULTS

- 250,200 were included in the analysis. Expressed variants and fusions are reported amongst nodules that were Afirma GSC-(S)uspicious or that had Bethesda V/VI cytology.
- 8,627 samples with TERT promoter mutation testing were analyzed.
- A subset of 17,000 GSC-S or B V/VI samples with available transcriptomic were analyzed for ERK and EMT signaling expression, the BRAF-RAS (BRS) score, and expression of the Human Molecular Signatures Database (MSigDB) hallmarks of cancer pathways.

TABLE 1

Proportion of Bethesda cytology categories by age group

- Patients ≥75 yrs had a higher rate of BIII cytology compared to those 20-54 yrs.
- Patients ≥75 yrs had lower rates of BV/VI cytology compared to those ≤74 yrs.

	Age groups				
	20-54	55-74	≥75		
Total	99,671	116,788	32,141		
B-III	75,632 (76%)	94,472 (80.9%)	26,044 (81%)*		
B-IV	16,847 (16.9%)	18,155 (15.5%)	5,136 (16%)		
B-V	3,416 (3.4%)	2,182 (1.9%)	535 (1.7%)**		
B-VI	3,779 (3.8%)	1,978 (1.7%)	426 (1.3%)**		

*p <0.01 compared to 20-54 yrs; **p <0.01 compared to ≤74 yrs

 Patients ≥75 years had larger median nodule size (2.4 cm [1.7-3.3]) compared to 55-74 years (2.2 cm [1.6-3.1]) and <54 years (2.2 cm [1.5-3.3]) (p <0.01).

TABLE 2 Afirma GSC-B call rates by age group (a) and sex (b)

(a)		Age groups	
	20-54	55-74	≥75
Total	92,479	112,627	31,180
GSC-B	56,015 (60.6%)	83,390 (74%)	23,648 (76%)*

*p <0.01 in ≥75 years compared to younger ages.

(b)		Age groups (Male)	
	20-54	55-74	≥75
Total	15,398	28,988	9,094
GSC-B	8,832 (57.4%)	19,584 (67.5%)	6,237 (68.6%)

	Age groups (Female**)		
	20-54	55-74	≥75
Total	77,058	83,613	22,084
GSC-B	47,169 (61.2%)	63,789 (76.3%)	17,409 (79%)

**p <0.01 in females compared to males.

- Patients ≥75 years had a higher Afirma GSC-(B)enign rate compared to younger age groups (Table 2a).
- Females had a higher rate of Afirma GSC-B than males in all age groups (Table 2b).

TABLE 3

Proportion of expressed variants and fusions in Afirma GSC-S (a) and BV/VI nodules (b)

Molecular alterations were detected at different rates across age groups in indeterminate nodules (BIII/IV) that are Afirma GSC-S (Table 3a) and BV/VI nodules (Table 3b).

(a)	Age groups	20-54	55-74	≥75	P values	
	Afirma GSC-S	36,449	29,235	7,530	20-54 vs ≥75	≥75 vs 55-74
	Total variant alterations	14,861 (40.77%)	10,206 (34.91%)	2,515 (33.4%)	<0.0001	0.01
	BRAFp.K601E	462 (1.27%)	285 (0.97%)	55 (0.73%)	0.0001	ns
	BRAFp.V600E	3,172 (8.7%)	1,969 (6.74%)	384 (5.1%)	<0.0001	<0.0001
	HRAS	3,257 (8.94%)	2,513 (8.6%)	607 (8.06%)	0.01	ns
	KRAS	637 (1.75%)	537 (1.84%)	129 (1.71%)	ns	ns
TS	NRAS	5,790 (15.89%)	3,578 (12.24%)	895 (11.89%)	<0.0001	ns
VARIANTS	DICER1	798 (2.19%)	293 (1%)	60 (0.8%)	<0.0001	ns
A A	EIF1AX	106 (0.29%)	204 (0.7%)	97 (1.29%)	<0.0001	<0.0001
	SPOP	193 (0.53%)	335 (1.14%)	88 (1.17%)	<0.0001	ns
	JAK2	22 (0.06%)	90 (0.3%)	66 (0.87%)	<0.0001	<0.0001
	PIK3CA	27 (0.07%)	50 (0.17%)	20 (0.26%)	<0.0001	ns
	TP53	47 (0.13%)	60 (0.2%)	27 (0.36%)	<0.0001	0.02
	TSHR	399 (1.09%)	377 (1.3%)	102 (1.35%)	ns	ns
	Total fusion alterations	2,549 (6.99%)	1,282 (4.39%)	266 (3.53%)	<0.0001	0.001
SZ	RET	356 (0.98%)	127 (0.43%)	29 (0.39%)	<0.0001	ns
FUSIONS	NTRK3	541 (1.48%)	165 (0.56%)	31 (0.41%)	<0.0001	ns
F	NTRK1	56 (0.15%)	29 (0.1%)	8 (0.11%)	ns	ns
	ALK	114 (0.31%)	50 (0.17%)	9 (0.12%)	0.005	ns
_	Total	3,066	2,538	750		
TERT	TERT+	46 (1.5%)	138 (5.4%)	81 (10.8%)	<0.0001	<0.0001
•	TERT+BRAFV600E	5 (0.2%)	32 (1.3%)	17 (2.3%)	<0.0001	0.06

(b)	Age groups	20-54	55-74	≥75	P values	
	BV/VI nodules	7,195	4,160	961	20-54 vs ≥75	≥75 vs 55-74
	Total variant alterations	4,934 (68.58%)	2,736 (65.77%)	599 (62.33%)	0.0001	0.05
	BRAFp.K601E	4 (0.06%)	9 (0.22%)	2 (0.21%)	ns	nd
	BRAFp.V600E	4,561 (63.39%)	2,426 (58.32%)	509 (52.97%)	<0.0001	0.002
	HRAS	68 (0.95%)	47 (1.13%)	12 (1.25%)	ns	ns
	KRAS	39 (0.54%)	26 (0.63%)	12 (1.25%)	0.01	ns
ZI	NRAS	131 (1.82%)	87 (2.09%)	31 (3.23%)	0.005	0.05
ARIANTS	DICER1	53 (0.74%)	24 (0.58%)	4 (0.42%)	ns	ns
X	EIF1AX	7 (0.1%)	9 (0.22%)	8 (0.83%)	<0.0001	0.007
	SPOP	21 (0.29%)	27 (0.65%)	10 (1%)	0.001	ns
	JAK2	3 (0.04%)	4 (0.1%)	2 (0.2%)	ns	ns
	PIK3CA	19 (0.26%)	28 (0.67%)	13 (1.35%)	<0.0001	ns
	TP53	10 (0.14%)	13 (0.31%)	6 (0.6%)	0.005	ns
	TSHR	39 (0.54%)	58 (1.39%)	8 (0.83%)	ns	ns
	Total fusion alterations	716 (9.95%)	241 (5.79%)	39 (4.06%)	<0.0001	0.04
SN	RET	328 (4.56%)	84 (2.02%)	8 (0.83%)	<0.0001	0.02
FUSIONS	NTRK3	151 (2.1%)	33 (0.79%)	3 (0.31%)	<0.0001	ns
3	NTRK1	35 (0.49%)	7 (0.17%)	1 (0.1%)	ns	ns
	ALK	29 (0.4%)	7 (0.17%)	1 (0.1%)	ns	ns
_	Total	1,275	723	175		
TERT	TERT+	45 (3.6%)	107 (14.8%)	55 (31.4%)	<0.0001	<0.0001
	TERT+BRAFV600E	38 (3%)	81 (11.2%)	41 (23.4%)	<0.0001	<0.0001

FIGURE 1

The activity expression of ERK and EMT and the BRS score by age groups

ERK and EMT activity expression is significantly lower with advancing age groups while there is no change in the BRS score (Figure 1).





FIGURE 2

MSigDB Hallmarks of Cancer Pathways in patients ≥75 years vs younger patients

Hallmarks where the confidence intervals does not cross the reference line are statistically significant (p < 0.05)

The expression of MSigDB hallmarks of cancer pathways are enriched in markers of UV response, KRAS signaling, MTORC1 signaling and immune signaling, amongst others, in those ≥75 years vs younger patients. Hedgehog, TGFβ, and EMT signaling are lower (Figure 2).



CONCLUSION

- Molecularly tested nodules from patients ≥75 years have an overall higher benign call rate and generally favorable variant/fusion profiles.
- In BV/VI nodules, there are higher levels of alterations suggestive of aggressive thyroid cancer (e.g. *TERT*p).
- Differential gene expression analysis suggests nodules in patients ≥75 years have increased inflammatory changes, while younger patients demonstrate more neoplastic pathways.
- Future studies should be done with histopathology detail and longer-term oncologic outcomes to see if improved prognostic features from the above-described expression analysis can personalize treatment for thyroid cancer in older patients.
- These data suggest one can increase confidence in more conservative management in older patients and may decrease the need for unnecessary surgery in this population.

REFERENCES

Wang Z, Vyas CM, Van Benschoten O, Nehs MA, Moore FD Jr, Marqusee E, Krane JF, Kim MI, Heller HT, Gawande AA, Frates MC, Doubilet PM, Doherty GM, Cho NL, Cibas ES, Benson CB, Barletta JA, Zavacki AM, Larsen PR, Alexander EK, Angell TE. Quantitative Analysis of the Benefits and Risk of Thyroid Nodule Evaluation in Patients ≥70 Years Old. Thyroid. 2018 Apr;28(4):465-471. doi: 10.1089/thy.2017.0655. Epub 2018 Apr 2. PMID: 29608439.
 Patel KN, Angell TE, Babiarz J, Barth NM, Blevins T, Duh QY, Ghossein RA, Harrell RM, Huang J, Kennedy GC, Kim SY, Kloos RT, LiVolsi VA, Randolph GW, Sadow PM, Shanik MH, Sosa JA, Traweek ST, Walsh PS, Whitney D, Yeh MW, Ladenson PW. Performance of a Genomic Sequencing Classifier for the Preoperative Diagnosis of Cytologically Indeterminate Thyroid Nodules. JAMA Surg. 2018 Sep 1;153(9):817-824. doi: 10.1001/jamasurg.2018.1153. PMID: 29799911; PMCID: PMC6583881.

© 2025 Veracyte, Inc. All rights reserved. AG1967.1.2508