

FNAprotect

Version 1.5 Revision Date 03/15/2016 Print Date 03/23/2016

SECTION 1. PRODUCT AND COMPANY IDENTIFICATION

Product name : FNAprotect

Manufacturer or supplier's details

Company : QIAGEN GmbH

QIAGEN Str. 1 D-40724 Hilden

Telephone : +49-02103-29-0

Responsible Department : QIAGEN Technical Service, QIAGEN Inc.,

27220 Turnberry Lane, Valencia, CA 91355, USA,

Tel.: 1-800-DNA-PREP (1-800-362-7737), Email: techservice-us@qiagen.com

E-mail : cpc@qiagen.com

addressResponsible/issuing

person

Emergency telephone : CHEMTREC

USA & Canada 1-800-424-9300

Recommended use of the chemical and restrictions on use

Recommended use : Laboratory chemicals

SECTION 2. HAZARDS IDENTIFICATION

GHS Classification

Skin irritation : Category 2
Serious eye damage : Category 1
Skin sensitization : Category 1
Acute aquatic toxicity : Category 2
Chronic aquatic toxicity : Category 2

GHS Label element

Hazard pictograms





Signal Word : Danger

Hazard Statements : H315 Causes skin irritation.

H317 May cause an allergic skin reaction.



FNAprotect

Version 1.5 Revision Date 03/15/2016 Print Date 03/23/2016

H318 Causes serious eye damage.

H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements : Prevention:

P280 Wear protective gloves/ protective clothing/ eye

protection/ face protection.

Response:

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and

easy to do. Continue rinsing.

P308 + P311 IF exposed or concerned: Call a POISON

CENTER or doctor/ physician.

Potential Health Effects

Aggravated Medical

Condition

Symptoms of Overexposure

: None known.

: No information available.

Carcinogenicity:

IARC No ingredient of this product present at levels greater than or

equal to 0.1% is identified as probable, possible or confirmed

human carcinogen by IARC.

ACGIH No ingredient of this product present at levels greater than or

equal to 0.1% is identified as a carcinogen or potential

carcinogen by ACGIH.

OSHA No ingredient of this product present at levels greater than or

equal to 0.1% is identified as a carcinogen or potential

carcinogen by OSHA.

NTP No ingredient of this product present at levels greater than or

equal to 0.1% is identified as a known or anticipated carcinogen

by NTP.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Hazardous ingredients

Chemical Name	CAS-No.	Concentration [%]
tetradecyltrimethylammonium oxalate	154858-16-9	>= 1 - < 5

SECTION 4. FIRST AID MEASURES

General advice : Move out of dangerous area.

Show this material safety data sheet to the doctor in

attendance.



FNAprotect

Version 1.5 Revision Date 03/15/2016 Print Date 03/23/2016

If inhaled : If unconscious place in recovery position and seek medical

advice.

If symptoms persist, call a physician.

In case of skin contact : Wash off immediately with soap and plenty of water while

removing all contaminated clothes and shoes.

If symptoms persist, call a physician.

In case of eye contact : Small amounts splashed into eyes can cause irreversible

tissue damage and blindness.

In the case of contact with eyes, rinse immediately with plenty

of water and seek medical advice.

Remove contact lenses. Protect unharmed eye.

If swallowed : If accidentally swallowed obtain immediate medical attention.

Rinse mouth with water.

Never give anything by mouth to an unconscious person.

Most important symptoms and effects, both acute and

delayed

: Causes skin irritation.

May cause an allergic skin reaction.

Causes serious eye damage.

No information available.

Notes to physician : No information available.

SECTION 5. FIRE-FIGHTING MEASURES

Suitable extinguishing media : Use water spray, alcohol-resistant foam, dry chemical or

carbon dioxide.

Specific hazards during fire

fighting

: Do not allow run-off from fire fighting to enter drains or water

courses.

Exposure to decomposition products may be a hazard to

health.

Hazardous combustion

products

Carbon monoxide, carbon dioxide and unburned

hydrocarbons (smoke).

Specific extinguishing

methods

: In the event of fire and/or explosion do not breathe fumes.

Special protective equipment

for fire-fighters

Wear self-contained breathing apparatus for firefighting if

necessary.

SECTION 6. ACCIDENTAL RELEASE MEASURES



FNAprotect

Version 1.5 Revision Date 03/15/2016 Print Date 03/23/2016

Personal precautions, protective equipment and emergency procedures Environmental precautions : Use personal protective equipment.

Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray.

: Prevent product from entering drains.

Prevent further leakage or spillage if safe to do so.

Methods and materials for containment and cleaning up

: Soak up with inert absorbent material (e.g. sand, silica gel,

acid binder, universal binder, sawdust).

Keep in suitable, closed containers for disposal.

SECTION 7. HANDLING AND STORAGE

Advice on safe handling : Do not breathe vapors/dust.

Avoid exposure - obtain special instructions before use.

Avoid contact with skin and eyes. For personal protection see section 8.

Smoking, eating and drinking should be prohibited in the

application area.

Dispose of rinse water in accordance with local and national

regulations.

Persons susceptible to skin sensitization problems or asthma, allergies, chronic or recurrent respiratory disease should not be employed in any process in which this mixture is being

used.

Conditions for safe storage : Keep container tightly closed in a dry and well-ventilated

place.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with workplace control parameters

Contains no substances with occupational exposure limit values.

Personal protective equipment

Hand protection

Remarks : The choice of an appropriate glove does not only depend on

its material but also on other quality features and is different from one producer to the other. Take note of the information given by the producer concerning permeability and break through times, and of special workplace conditions

(mechanical strain, duration of contact).

Material : Protective gloves

Eye protection : Tightly fitting safety goggles

Wear face-shield and protective suit for abnormal processing

problems.

Do not wear contact lenses.

Ensure that eyewash stations and safety showers are close to



FNAprotect

Version 1.5 Revision Date 03/15/2016 Print Date 03/23/2016

the workstation location.

Skin and body protection : Choose body protection according to the amount and

concentration of the dangerous substance at the work place.

acid-resistant protective clothing
Footwear protecting against chemicals

Hygiene measures : Keep away from food and drink.

Wash hands before breaks and at the end of workday. Ensure adequate ventilation, especially in confined areas.

Keep working clothes separately.

Avoid contact with the skin and the eyes. When using do not eat, drink or smoke.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Color : No data available

Odor : characteristic

Odor Threshold : No data available

pH : No data available

Melting point/range : No data available

Boiling point/boiling range : No data available

Evaporation rate : No data available

Burning rate : No data available

Upper explosion limit : No data available

Lower explosion limit : No data available

Vapor pressure : No data available

Relative vapor density : No data available

Relative density : No data available

Density : No data available

Solubility(ies)

Flash point

Water solubility : No data available

Solubility in other solvents : No data available

Partition coefficient: n-

octanol/water

: No data available

: No data available

Autoignition temperature : not determined

Thermal decomposition : No data available

Viscosity



FNAprotect

Version 1.5 Revision Date 03/15/2016 Print Date 03/23/2016

Viscosity, dynamic : No data available
Viscosity, kinematic : No data available
Explosive properties : No data available
Oxidizing properties : No data available

SECTION 10. STABILITY AND REACTIVITY

Reactivity : No decomposition if stored and applied as directed.

Chemical stability : No decomposition if stored and applied as directed.

Possibility of hazardous

reactions

: Stable under recommended storage conditions.

Hazardous decomposition products formed under fire

conditions.

Keep away from oxidizing agents, and acidic or alkaline

products.

Conditions to avoid : No data available

Incompatible materials : No data available

Hazardous decomposition

products

: No decomposition if stored and applied as directed.

SECTION 11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : No data available

Acute inhalation toxicity : No data available

Acute dermal toxicity : No data available

Skin corrosion/irritation

Product:

May cause skin irritation and/or dermatitis.

Serious eye damage/eye irritation

Product:

Remarks:



FNAprotect

 Version 1.5
 Revision Date 03/15/2016
 Print Date 03/23/2016

May cause irreversible eye damage.

Respiratory or skin sensitization

Product:

Remarks:

Causes sensitization.

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT-single exposure

No data available

STOT-repeated exposure

No data available

Aspiration toxicity

No data available

Further information

No data available

SECTION 12. ECOLOGICAL INFORMATION

Ecotoxicity

Product:

Toxicity to fish : No data available

Toxicity to algae : No data available

Toxicity to bacteria : No data available

Persistence and degradability

No data available

Bioaccumulative potential

Product:



FNAprotect

Version 1.5 Revision Date 03/15/2016 Print Date 03/23/2016

Bioaccumulation : No data available

Partition coefficient: n-

octanol/water

: No data available

Mobility in soil

No data available

Other adverse effects

No data available

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : The product should not be allowed to enter drains, water

courses or the soil.

Send to a licensed waste management company.

Dispose of as hazardous waste in compliance with local and

national regulations.

Contaminated packaging : Dispose of as unused product.

Do not re-use empty containers.

SECTION 14. TRANSPORT INFORMATION

IATA-DGR

UN/ID No. : 3082

Proper shipping name : Environmentally hazardous substance, liquid, n.o.s.

(TETRADECYLTRIMETHYLAMMONIUM OXALATE)

Class : 9
Packing group : III
Labels : 9

IMDG-Code

UN number : 3082

Proper shipping name : ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID,

N.O.S.

(TETRADECYLTRIMETHYLAMMONIUM OXALATE)

Class : 9
Packing group : III
Labels : 9
EmS Code : ,
Marine pollutant : no

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

No data available



FNAprotect

Version 1.5 Revision Date 03/15/2016 Print Date 03/23/2016

Domestic regulation

49 CFR

UN/ID/NA number : 3082

Proper shipping name : Environmentally hazardous substance, liquid, n.o.s.

(TETRADECYLTRIMETHYLAMMONIUM OXALATE)

Class : 9
Packing group : III
Labels : 9
ERG Code : 171

Marine pollutant : yes (TETRADECYLTRIMETHYLAMMONIUM OXALATE)

Special precautions for user:

For personal protection see section 8.

SECTION 15. REGULATORY INFORMATION

OSHA Hazards : Skin sensitizer, Corrosive to skin, Corrosive to eyes, Corrosive

to respiratory system.

EPCRA - Emergency Planning and Community Right-to-Know

CERCLA Reportable Quantity

This material does not contain any components with a CERCLA RQ.

SARA 304 Extremely Hazardous Substances Reportable Quantity

This material does not contain any components with a section 304 EHS RQ.

SARA 311/312 Hazards : Acute Health Hazard

SARA 302 : No chemicals in this material are subject to the reporting

requirements of SARA Title III, Section 302.

SARA 313 : This material does not contain any chemical components with

known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

US State Regulations

California Prop 65 This product does not contain any chemicals known to the State

of California to cause cancer, birth, or any other reproductive

defects.

SECTION 16. OTHER INFORMATION

The information provided in this Material Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is



Print Date 03/23/2016

FNAprotect

Version 1.5 Revision Date 03/15/2016

designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.